

**A.R.C.  
LASER**

# C-LAS ARM

Standard in surgery



**30 W**

**IDEAL IN  
HOSPITAL**



Laser handpiece, 6 inch tip



C-PRO  
Micromanipulator



**Focus  
200 to  
500 mm  
up to  
100 µm**

## C-LAS ARTICULATED ARM, 30 WATT

Amount	Order #	Description
1 x	LS02100-arm	CO <sub>2</sub> Laser 10.6 µm



## Use in OT (Operation Theatre)

### ACCESSORIES

PCS=1 | HS11024



CO<sub>2</sub> Laser handpiece: 6" tip

PCS=1 | ZU12000



C-Pro (micromanipulator)


PCS=1 | ZU11029-b



Medical Cart - Solo

## EYE SAFETY

AS01023 safety goggle PCS=1



## ViroVac Smoke Evacuation ZU01049



### ACCESSORIES ViroVac

PCS=10 | ZU01048



Flexible tubing 3.1m, d=22mm

PCS=10 | ZU01058s



Flexible tube, l=60cm, sterile

PCS=2 | ZU01050



Filter ViroSafe (spare part)

PCS=1 | KB01329



Remote cable ViroVac