

HNO

Endoskopie



BlueMed AG
Martinsbruggstrasse 65
9016 St.Gallen

T 044 741 52 32
info@blue-med.ch

www.blue-med.ch
shop.blue-med.ch





ENT (Ear, Nose, Throat)

otoscope | HD | autoclavable






diameter / length	0° 		30° 			
2.7 / 50	10-0032-00		10-0033-00			
3.0 / 110	10-0653-00		10-0654-00			
3.0 / 141	10-0696-00		10-0697-00			
4.0 / 50	10-0034-00		10-0035-00			

sinuscope semiflex | wide angle HD | autoclavable | 4.0 mm is semiflex generation 2.0



semiflex

diameter / length	0° 		30° 	45° 	70° 	
2.7 / 110	10-0036-00SF		10-0037-00SF			
2.7 / 175	10-0039-00SF		10-0040-00SF			
4.0 / 175	10-0045-00SF		10-0046-00SF	10-0086-00SF	10-0047-00SF	Semiflex generation 2.0



sinuscope | wide angle HD | autoclavable

diameter / length	0° 		30° 	45° 	70° 	90° 
1.9 / 69	10-0422-00		10-0416-00			
2.4 / 98	10-0264-00		10-0265-00			
2.7 / 110	10-0036-00		10-0037-00	10-0418-00	10-0038-00	
2.7 / 141	10-0466-00		10-0467-00	10-0486-00		
2.7 / 175	10-0039-00		10-0040-00	10-0387-00	10-0041-00	
2.7 / 187.5	10-0042-00		10-0043-00		10-0044-00	
3.0 / 110	10-0573-00		10-0574-00			
3.0 / 141	10-0576-00		10-0577-00			
4.0 / 140	10-0150-00		10-0149-00		10-0151-00	
4.0 / 175				10-0086-00		
4K 4.0 / 175 (4K)	10-2045-00		10-2046-00		10-2047-00	



microlaryngoscope | HD | autoclavable

diameter / length				70° 	90° 
4.0 / 175				10-0090-00	10-0475-00
5.8 / 188				10-0219-00	



laryngoscope | HD | lateral illumination | autoclavable

diameter / length				70° 	90° 
8.0 / 190				10-0177-00	
10.0 / 190				10-0476-00	10-0477-00
flat 12.4x7.2 / 184				10-0626-00	10-0635-00


telemetryngoscope | New Focus Design HD | lateral illumination; extension: 5x / 7x | soakable

diameter / length				70° 	90° 
10.0 / 200				10-0468-00	10-0442-00
flat 12.4x7.2 / 200				10-0768-00	10-0742-00

bronchoscope | HD | autoclavable

diameter / length	0° 		30° 			
2.9 / 303	10-0618-00		10-0619-00			
2.9 / 365	10-0383-00		10-0384-00			
5.5 / 490	10-0050-00		10-0051-00			

oesophagoscope | HD | autoclavable

diameter / length	0° 					
5.5 / 530	10-0049-00					

For any further requirements, please contact us. We are able to design any endoscopic product specific to your needs. Semiflex is protected as utility model and patented.

TV Adapter

HD design | optimised for HD-camera | c-mount
| catch asap standard | with special focus ring for HD camera head

focus	soakable
f 14	10-0480-00
f 21	10-0488-00
f 25	10-0445-00
f 28	10-0448-00

vario-zoom | parfocal | c-mount | catch asap standard

focus	soakable
f 15 – f 35	10-0238-00

zoom | c-mount | catch asap standard

focus	soakable
f 15 – f 25	10-0060-00

standard design | c-mount | catch asap standard

focus	soakable
f 15	10-0055-00
f 21	10-0052-00
f 25	10-0110-00
f 28	10-0053-00
f 35	10-0054-00
f 40	10-0074-00

short design | c-mount | catch asap standard | height: 38.5mm

focus	soakable
f 15	10-0155-00
f 21	10-0152-00
f 28	10-0153-00

beamsplitter | splitting 90/10 | c-mount | catch asap standard

focus	soakable
f 18	10-0132-00
f 25	10-0131-00
f 30	10-0130-00

angled 90° | c-mount | catch asap standard

focus	soakable
f 18	10-0146-00
f 25	10-0137-00
f 30	10-0136-00



For any further requirements, please contact us. We are able to design any endoscopic product specific to your needs.



10'000+
verfügbare Produkte



BlueMed AG
Martinsbruggstrasse 65
9016 St.Gallen

T 044 741 52 32
info@blue-med.ch

www.blue-med.ch
shop.blue-med.ch